

YOUR RECEIVER MAINTENANCE KIT.

ALL YOU NEED TO MANAGE YOUR HEARING HEALTHCARE.

signia

Life sounds brilliant.



The Receiver Maintenance Kit

ESSENTIAL FOR YOUR HEARING HEALTH

Styletto has transformed the look of hearing aids with its new slim, stylish, unique hearing aid design which is attracting new and existing customers alike.

As with all hearing aids, Styletto requires regular maintenance and care which not only maintains the maximum performance of the hearing aid but also ensures the preservation of your customers hearing health.

To help you manage your Styletto customers hearing healthcare, Signia provides a Receiver Maintenance Kit mailer box which you can supply to your customers on a regular basis along with with the required domes/sleeves and miniReceiver wax guards.

We recommend exchanging Click Domes/Sleeves every 1-2 months or sooner if there are cracks or other changes. Replace wax guards every 1-2 months or sooner if there is noticeable wax build up.


For more information and to order your Receiver Maintenance Kit mailer boxes, please contact your Account Manager.

How it works?

- Sivantos provides you with the Receiver Maintenance Kit mailer boxes - ask your Account Manager to order your Receiver Maintenance Kit mailer boxes
- You order the required click domes/sleeves and miniReceiver wax guards for your client(s)
- You send the Domes/Sleeves and wax guards to your client(s) using the Receiver Maintenance Kit mailer box and charge accordingly
- We recommend you send your clients their Receiver Maintenance Kit every 2-3 months

Part numbers for ordering

Click Domes	Type	Units/Box	Price
Open	4 mm Open	6	£7.25
	8 mm Open	6	£7.25
	10 mm Open	6	£7.25
Semi-open	Semi-Open	6	£7.25
	6 mm Closed	6	£7.25
Closed	8 mm Closed	6	£7.25
	10 mm Closed	6	£7.25
	8/10 mm Double	6	£7.25
Double	10/12 mm Double	6	£7.25

Click Sleeves	Size	Type	Price
	XS	Closed	£8.70
	S	Closed	£8.70
	M	Closed	£8.70
	L	Closed	£8.70

Accessories	Description	Price
	miniReceiver Wax Guard (10 pcs)	£3.50

Other products for ongoing receiver care



Click Dome Removal Tool
£3.00



PerfectDry Lux
£50.00



Perfect Clean
£90.00



Hearing Aid Cleaning Set
£9.00

Set Includes Bag, 1x Tool Set, 1x Drying Cup, 1x Cleaning Spray, 4x Drying Capsules, 5x Cleaning Wipes, 1x Cleaning Cloth

signia-pro.co.uk/styletto