

Insio

signia

Life sounds brilliant.

REPLICATING NATURE.

The new Insio hearing aids from Signia.
The Bluetooth enabled in-the-ear hearing aids
for best speech comprehension.



signia-hearing.com/insio

Full control in your hands.



myControl App



StreamLine TV

StreamLine Mic



miniPocket



Telephone

Stream phone calls directly from your iPhone (and with **StreamLine Mic** from all smartphones) to your hearing aids for complete convenience.



Music

Listen to your favorite music streamed directly into your hearing aids from your Apple device (and with **StreamLine Mic** from any device with Bluetooth) for an outstanding sound experience.



TV

Stream audio from your television directly into your hearing aids with **StreamLine TV**.



Remote control

Change Insignia's settings remotely using the **myControl App™***.

Find out more about the benefits and convenience of connectivity:
www.signia-hearing.com/connectivity

Experience the joy of natural hearing.



Speech in noise

Thanks to its **Speech Quality**, Signia's Insio enables you to hear speech clearly in all situations through sophisticated microphones and advanced noise suppression technologies. This makes it easier for you to understand your conversation partner, even in very noisy environments such as a restaurant or a conference.



Conversation while moving

The **3D Classifier** of Insio provides optimal hearing of your own voice and your acoustic environment even when moving - it automatically adjusts the settings based on the **myControl App's motion data**.



Crystal clear sound

Insio's **Sound Clarity** offers the most natural, lifelike sound in all environments. This lets you enjoy crystal clear sound, even in the loudest environments.

Your partner for confidence and natural sound.



Most natural hearing experience

With Insignia by Signia, you can understand your conversation partner even in difficult hearing environments while it offers you perfectly balanced sound.



Dedicated push button

Via a push button, you can change your hearing programs or the volume manually.



For a broad range of hearing loss


Mild, moderate or severe hearing loss, Insignia tackles all these challenges and gives you back the natural sound of life.



The new Insio.

Life is better spent together. Signia's latest Insio™ hearing aids are tailor-made to suit your ear anatomy and hearing loss so you can enjoy spending time with those closest to you without having to worry about your hearing and appearance.

Through its incredibly discreet design and exceptional sound quality, the new Insio gives you the self-confidence to savor the important things in life. Plus, Bluetooth connectivity in selected models provides carefree direct audio streaming of TV, music and phone calls.




Crystal clear sound with a natural flow.

From animatedly discussing music after a concert to enjoying the thrills of a new city, Signia Insignia lets you delight in the world of sound.

Conversations require concentration and effort wherever ambient sounds can distract us from hearing what matters. Insignia is expertly equipped to tackle these challenges. It replicates the natural way we localize sound, focusing on your conversation partner while filtering unwanted background noise. You hear most clearly what you want to hear. All day long. In every situation.

Please note:

Features may vary depending on the configuration of the individual hearing aid. Your hearing care professional will be happy to assist you.




Enjoy direct streaming of phone calls, music, and TV.

Stream music in high quality stereo straight into your ears, discreetly adjust the volume to suit your wishes while watching your favorite TV show, or hear a phone call at the push of a button.

With Insio you have discreet and convenient remote control and premium Bluetooth connectivity for audio streaming from various sources at your fingertips - all you need is a smartphone.

More about Insio.



Further facts and figures.



High-quality audio streaming


Thanks to their Bluetooth connectivity, Signia's Insio ITC and ITE models let you fully enjoy convenient and carefree direct streaming of TV, music and phone calls.



Signia tinnitus solutions

Tune out tinnitus with the unique built-in tinnitus therapies from Signia. Ask your hearing care professional which solution is best for you.

Actual size in cm



Available colors:



Beige



Tan



Brown



Dark
Brown




Mocha




myHearing App.

Your myHearing App



Your hearing care professional



Personal care, anywhere.

Your relaxed journey back to effortless hearing is even easier with the myHearing App* for your smartphone. It lets you familiarize yourself with your new hearing aids in as little time as possible. You can follow special listening exercises, consult the user guide and connect to your hearing care professional on the go.



Download on the
App Store



GET IT ON
Google Play

* Compatible with Android and iOS. The apps are free of charge.

Find out more on myhearing.com or ask your hearing care professional for more information on Insignia hearing aids, useful accessories and the myHearing App.

The information in this document contains general descriptions of the technical options available, which do not always have to be present in individual cases and are subject to change without prior notice

Android and Google Play are trademarks of Google Inc.
Apple App Store is a trademark of Apple Inc.

iPhone is a trademark of Apple Inc., registered in the U.S. and other countries.

Legal Manufacturer

Signia GmbH
Henri-Dunant-Strasse 100
91058 Erlangen
Germany

Signia GmbH is a Trademark Licensee of Siemens AG.

Order No. 03637-99C1-7600 | Printed in Germany
© 08.2018 Signia GmbH

Find us on:

signia-hearing.com



facebook.com/signiahearing



youtube.com/signiahearing